



D 4000

D 4000 is a complete digital system for acquisition and processing of digital images obtained from radiography and fluoroscopy studies

Main applications

- Standard dynamic procedures:
 - G.I. studies
 - Uro-genital
 - Artherio-phlebography
- Angiographic procedures:
 - Cerebral
 - Pulmonary
 - Lower-upper extremity
 - Coronary
 - Kidney
 - Liver

Main features

- 12" and 16" Image intensifiers
- 1024x1024 CCD digital camera
- 12 bit - 4096 gray levels contrast resolution
- Dynamic acquisition up to 30 frames/sec
- Large series of image dedicated processing features
- Optional angiographic module

Technical data

Image intensifier	12" , 3 fields
	16" , 4 fields
Optical system	The compact type optical system is specifically designed to transfer the output image from image intensifier to CCD camera. It is complete of remote controlled iris diaphragm and photodiode optical sensor for dose control.
TV system	The TV camera uses a large area CCD sensor: 13.1x13.1 mm. Resolution:1024x1024x12 bit (4096 grey levels) Dynamic range: > 66db Progressive scan rate up to 30 frames /s. in pulsed fluoroscopy and from single up to 15 frames/s in fluorography
	CCD sensor offer significant advantages compared to conventional pick-up tube: - resolution limited only by pixel number and not further degraded by the electron optics; - uniform resolution over the entire image without the typical reduction at the border; - low readout noise that allow the same signal-to-noise ratio is achieved with a lower signal, i.e. at a lower dose.
18" LCD Monitor	1280x1024 pixel
PC Hardware	The system is set up by a main controller Motorola 68332, a real time image processor and a host computer with CPU Intel. The computer is complete of 512 MB RAM memory for image acquisition and 36 GB hard disc for images storing. The operative system is Windows® NT 4.0 assuring easy and friendly user interface.

Acquisition modes and real time processing

Continuous fluoroscopy	<ul style="list-style-type: none">- matrix 1024x1024,- visualization at 30 frames/s,- image saving on hard by means a proper command,- last image hold,- noise reduction with motion sensitivity,- automatic windowing,- multistep edge enhancement,- horizontal and vertical digital image reversal.
Pulsed digital fluoroscopy	<ul style="list-style-type: none">- matrix 1024x1024,- acquisition rate up to 30 frames/s,- automatic image saving on hard disc,- last image hold,- automatic windowing,- multistep edge enhancement,- horizontal and vertical digital image reversal.
Digital Radiography	<ul style="list-style-type: none">- matrix 1024x1024,- acquisition rates single frame and from 0.5 up to 15 frames/s,- tomography,- direct image saving on hard disc,- automatic windowing,

- multistep edge enhancement,
- horizontal and vertical digital image reversal.

Processing

Contrast, brightness, gamma control
 Spatial filters
 Enlargements
 Overlay and measurements
 Electronic diaphragms
 Overview
 Harmonization

Image display

18" LCD High resolution monochromatic monitors with high brightness, contrast and resolution
 For main or reference application in the diagnostic room 1 or 2 additional monitors are possible on request

Post processing

- Equalization histogram in the interest area;
- High speed spatial filters (Smooth, Sharp) with variable Kernel size (from 3x3 to 11x11);
- Overview of the interest images: 4, 9, 16, 1+5, 1+7;
- Grey scale inversion;
- Window, level and gamma correction;
- Multistep edge enhancement
- Rectangular and circular automatic electronic collimators;
- Zoom with magnification factor from 1.2 to 3;
- Magnification lens, magnification factor from 1.25 to 3;
- Horizontal and vertical image shift;
- Left/right and up/down image flip;
- Cine loop;
- Framesage overlay of pre-programmed or user

Memory

RAM 512 MB

Hard Disk:

36 GB images

System Options

Additional Monitor 18"
 Harddisc 72 GB
 Kit Angiography
 Infrared remote control
 Software for „Image Stitching“
 DICOM Basis Interface- Image storage
 DICOM Print
 DICOM worklist
 DICOM CD ROM drive
 DICOM package: Storage, Print, Worklist
 Laser printer Digital Interface
 Scan down Konverter SVGA/PAL
 Additional console

External connections

RS 422 Laser printer Digital Interface
 Ethernet TCP/IP Network connectivity with DICOM 3 protocol