


RADRAY RS 10 /20

Wall stand fitting to FLUORORAY T systems



RADRAY RS 10 and RS 20 are Bucky wall stands matching the design line of FLUORORAY T fluoroscopy systems

Consisting of:

- Column stand for wall and floor mount
- Bucky with movable grid for cassette sizes up to 35x43 cm, counterbalanced with horizontal movement
- Electro-mechanical brakes for all movements
- Complete with grid and prepared to accommodate the AEC chamber
- RS 10 standard version / left loading
- RS 20 additional tilting of Bucky

| Technical data | |
|---|-------------------------------------|
| RADRAY RS 10 | |
| Vertical movement of Bucky | 136,5 cm, counter balanced |
| Brakes | Mechanical brakes |
| Minimum floor- center-line distance | Appr. 46 cm |
| Weight | 170 Kg |
| RADRAY RS 20 | |
| Tilting of Bucky | -20 to 90 degrees |
| Rotation of Bucky | 90 degrees |
| Minimum floor- center-line distance | Appr. 46 cm |
| Minimum floor- Bucky distance in horizontal position of the Bucky | Appr. 68,5 cm |
| Bucky | |
| Size of cassettes | 13x18 to 35x43 cm |
| Grid | 90 Line pair/inch R: 12:1; F=180 cm |
| Options and accessories | |
| Automatic cassette size sensing | |
| Wall-floor mounting frame | |
| Floor mount-only- | |
| AEC chamber | |
| Patient handle | |

CE 0473

STARMED International Medical Systems GmbH • An der Ronne 7 • D-50259 Pulheim
 Fon: +49 (0) 2234 658 208 • Fax: +49 (0) 2234 658 210 • E-Mail: info@starmedpro.com • http://www.starmedpro.com